NOTES:

SUBSTRATE:


GRADE A FINE ANNEALED SCHOTT: N-SF11 785/258

- 2. ROHS COMPLIANT
- 3. CENTERING TOLERANCE (AT 587.6nm): BEAM DEVIATION (HALF ANGLE): <1 ARCMIN
- 4. COATING (APPLY ACROSS COATING APERTURE)

\$1 & \$2: NIR II $R(ABS) \le 1.5\%$ FROM 750-800nm @ 0° AOI $R(ABS) \le 1.0\%$ FROM 800-1550nm @ 0° AOI $R(AVG) \le 0.7\%$ FROM 750-1550nm @ 0° AOI



- 6. POWER, IRREGULARITY, AND SURFACE QUALITY SPECIFICATIONS APPLY ACROSS CLEAR APERTURE
- 7. FOCAL LENGTH (EFL): -15.00mm±1% BACK FOCAL LENGTH (BFL): -16.68mm
- 8. PROTECTIVE BEVEL AS NEEDED
- 9. DESIGN WAVELENGTH: 587.6nm



SECTION A-A

FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING

	\$1	\$2			
SHAPE	CONCAVE	PLANO			
RADIUS	11.77	INFINITY			
SURFACE QUALITY	40 - 20	40 - 20			
MIN CLEAR APERTURE	Ø11.00 Ø11.00				
MIN COATING APERTURE	Ø11.00	Ø11.00			
POWER AT 632.8nm	3.00 RINGS 3.00 RINGS				
IRREGULARITY AT 632.8nm	GULARITY AT 632.8nm 0.50 RINGS 0.50 F				

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

		Edmund Optics ®		
THIRD ANG PROJECTIO		TITLE	12mm Diameter x -15 FL, NIR II Coated, Plano-Concave Lens	
ALL DIMS IN	mm	DWG NO	67992	SHEET 1 OF 1