




**FOR INFORMATION ONLY:  
DO NOT MANUFACTURE  
PARTS TO THIS DRAWING**



SHORTPASS FILTER

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE  
DIMENSIONS ARE FOR REFERENCE ONLY

REV. A	VALUE
CUT-OFF WAVELENGTH	450±2nm
TRANSMISSION WAVELENGTH	250 - 430nm
REJECTION WAVELENGTH	470 - 545nm
TRANSMISSION	>85
SURFACE QUALITY	40-20
COATING TYPE	TRADITIONAL COATED



THIRD ANGLE PROJECTION		TITLE	M25.5 x 0.5, 450nm MOUNTED SHORTPASS FILTER
ALL DIMS IN	mm	DWG NO	66640
			SHEET 1 OF 1