






**FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING**



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE
DIMENSIONS ARE FOR REFERENCE ONLY

	S1	S2	 Edmund Optics®			
SHAPE	PLANO	PLANO				
SURFACE QUALITY	60-40	60-40				
REFLECTION WAVELENGTH	400 - 430 496 - 510 & 579 - 596	400 - 430 496 - 510 & 579 - 596				
TRANSMITTED WAVEFRONT RMS	1λ	1λ				
BEVEL	PROTECTED AS NEEDED	PROTECTED AS NEEDED				
			ALL DIMS IN	mm	DWG NO	87278
			TITLE			Ø25mm, 444 521 & 608nm, MULTI-EDGE DICHROIC FILTER
						SHEET 1 OF 1