


FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING







-	- 5.00±0.10

	VALUE			
CENTER WAVELENGTH (nm)	355nm			
FULL WIDTH HALF MAX (nm)	10nm			
MINIMUM TRANSMISSION	≥85			
BLOCKING	200 - 1200			
SURFACE QUALITY	80-50			
CONSTRUCTION	MOUNTED IN BLACK ANODIZED ALUMINUM RING			
COATING TYPE	HARD COATED			

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

				Edmund Optic	\$S®
_	THIRD ANGLE PROJECTION		TITLE	OD4, Ø25mm, 355nm x 10nm, HC, BP I	FILTER
_	ALL DIMS IN	mm	DWG NO	34491	SHEET 1 OF 1