

1. SUBSTRATE: N-SF11 2. COATING (APPLY ACROSS COATING APERTURE) S1: NONE S2: NONE /3. FINE GRIND SURFACE 4. RoHS COMPLIANT — 3.00±0.05 – 3.00±0.05 S2--S2 /3\-3 2.12 NOMINAL 45°0'±0°2' --S1

NFORMA FO DO NOT MANUFACTURE PARTS TO THIS DRAWING

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

REV. A	S1	S2					
SHAPE	PLANO	PLANO				Edmund Optic	۲ د ا
SURFACE QUALITY	20-10	20-10					<i>י</i> כי
CLEAR APERTURE	3.82 x 2.70	2.70 x 2.70			TITLE	3mm, Uncoated, N-SF11 Micro Right Angle Prism	
COATING APERTURE	N/A	N/A					
SURFACE FLATNESS	λ/4	λ/4					
BEVEL	PROTECTIVE AS NEEDED	PROTECTIVE AS NEEDED	ALL DIMS IN	mm	DWG NO	45525	SHEET 1 OF 1