NOTES:


- 1. SUBSTRATE MATERIAL: Fused Silica
- 2. SURFACE S2 TO BE PARALLEL WITH S1 TO 30 arcsec
- 3. COATING (APPLY ACROSS CLEAR APERTURE)


S1: UV Enhanced Aluminum R(avg) >89% @ 250 - 450nm R(avg) >85% @ 250 - 700nm

S2: NONE

4. FINE GRIND SURFACE

5. Rohs Compliant





FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

REV. A	\$1	\$2		TITLE	25mm Diameter UV Enhanced Aluminum Coated, λ/10 Mirror	
SURFACE QUALITY	20-10	COMMERCIAL POLISH	THIRD ANGLE PROJECTION			
SURFACE FLATNESS	λ/10	NA			Coalea, A, 10 Millor	CUEET
BEVEL	PROTECTIVE AS NEEDED	PROTECTIVE AS NEEDED	ALL DIMS IN mm	DWG NO	68315	SHEET 1 OF 1

