

NOTES:



1. SUBSTRATE MATERIAL:
Fused Silica
2. COATING (APPLY ACROSS CLEAR APERTURE)

S1: Enhanced Aluminum
R(avg) >95% @ 450 - 650nm

S2: NONE

3.

FINE GROUND SURFACE
4. PARALLELISM: 30 ARCSEC
5. RoHS COMPLIANT



**FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING**

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE
DIMENSIONS ARE FOR REFERENCE ONLY

REV. A	S1	S2	<div><div><div><div></div></div></div><div><div>Edmund Optics®</div></div></div>			
SHAPE	PLANO	PLANO	<div>TITLE</div> 30mm Square Enhanced Aluminum Coated, λ/10 Mirror			
SURFACE QUALITY	20-10	N/A				
SURFACE FLATNESS	λ/10	COMMERICAL POLISH	<div>DWG NO</div> 84412			
CLEAR APERTURE	90%	N/A				
BEVEL	PROTECTIVE AS NEEDED	PROTECTIVE AS NEEDED	<div>ALL DIMS IN</div> mm	<div>SHEET</div> 1 OF 1		