

FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING





SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

Edmund Optics®

EXPANSION POWER	5X	THIRD ANGLE PROJECTION		TITLE	5X, 532nm Scorpii™ Nd:YAG Beam Expander	
DESIGN WAVELENGTH DWL (nm)	532					ander
INPUT APERTURE DIAMETER(mm)	3.5		1			SHEET
EXIT APERTURE DIAMETER(mm)	10.0	ALL DIMS IN	mm	DWG NO	37056	1 OF 1