




**FOR INFORMATION ONLY:  
DO NOT MANUFACTURE  
PARTS TO THIS DRAWING**



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE  
DIMENSIONS ARE FOR REFERENCE ONLY

NUMERICAL APERTURE NA	0.75
MAGNIFICATION	20X
WORKING DISTANCE (mm)	0.6
COVER GLASS THICKNESS	0.17
FIELD LENGTH FL (mm)	9
FIELD NUMBER	26.5
RoHS	COMPLIANT

THIRD ANGLE PROJECTION 		 <b>Edmund Optics®</b>		
TITLE		Olympus UPLSAPO 20X Objective		
ALL DIMS IN	mm	DWG NO	36064	SHEET 1 OF 1