





**FOR INFORMATION ONLY:  
DO NOT MANUFACTURE  
PARTS TO THIS DRAWING**



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE  
DIMENSIONS ARE FOR REFERENCE ONLY

NUMERICAL APERTURE NA	0.95
MAGNIFICATION	40X
WORKING DISTANCE (mm)	0.18
COVER GLASS THICKNESS	0.11 - 0.23
FIELD LENGTH FL (mm)	4.5
FIELD NUMBER	26.5
RoHS	COMPLIANT

 <b>Edmund Optics®</b>			
THIRD ANGLE PROJECTION 		TITLE	Olympus UPLSAPO 40X2 Objective
ALL DIMS IN	mm	DWG NO	36065
			SHEET 1 OF 1