



**FOR INFORMATION ONLY:  
DO NOT MANUFACTURE  
PARTS TO THIS DRAWING**



EFFECTIVE FOCAL LENGTH EFL (mm)	15.37
MAGNIFICATION	10X
MOUNTING THREAD	0.7965-36 TPI, 55° Whitworth
NUMERICAL APERTURE NA	0.25
WORKING DISTANCE (mm)	1.5
LENGTH EXCLUDING THREADS (mm)	44
MAXIMUM DIAMETER (mm)	24
TYPE	Plan
RoHS	NOT COMPLIANT

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE  
DIMENSIONS ARE FOR REFERENCE ONLY

**EO**® **Edmund Optics**®



THIRD ANGLE  
PROJECTION

TITLE

10X DIN Plan Commercial Grade Objective

ALL DIMS IN

mm

DWG NO

43907

SHEET  
1 OF 1