

FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING





COMPATIBLE TUBE LENS	Focal Length: 200mm
MAGNIFICATION	10X
MOUNTING THREAD	M26 x 36 TPI
NUMERICAL APERTURE NA	0.28
WORKING DISTANCE (mm)	33.5
MOUNTING THREADS	M26 x 36 TPI
FOCAL LENGTH FL (mm)	20
TYPE	Plan APO
FIELD OF VIEW, 1" SENSOR	1.28 x 0.96mm
FIELD OF VIEW, 1/2" SENSOR	0.64 x 0.48mm
FIELD OF VIEW, 2/3" SENSOR	0.88 x 0.66mm
FIELD OF VIEW, 24 DIAMETER FIELD EYEPIECE (mm)	2.4
DEPTH OF FOCUS (µm)	3.6
RESOLVING POWER (µm)	1
RoHS	NOT COMPLIANT

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

