

NOTES:

- CONSTRUCTION:
 Aluminum Housing, Blue Tempered Spring Steel Leaves
- 2. NUMBER OF LEAVES: 20
- 3. ROHS: COMPLIANT





SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

Edmund Optics®

FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING

THIRD ANGLE PROJECTION		TITLE	180mm Outer Diameter, Iris Diaphragm	
ALL DIMS IN	mm	DWG NO	32608	SHEET 1 OF 1