


NOTES:

- COATING (APPLY ACROSS COATING APERTURE)
Protected Gold (700-11000nm)
R(AVG) >96% @ 700 - 11000nm




SECTION A-A

**FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING**

SPECIFICATIONS SUBJECT TO
CHANGE WITHOUT NOTICE
DIMENSIONS ARE FOR
REFERENCE ONLY

EO® Edmund Optics®

THIRD ANGLE PROJECTION 		TITLE	Infinite Conjugate, Gold Coated, 20X/0.33NA RefIX Objective	
ALL DIMS IN	mm	DWG NO	39143	SHEET 1 OF 1