## NOTES:


1. CONSTRUCTION:

Stainless Steel (400 Series) Top and Bottom Plates, Plated Steel Honeycomb Core, High Pressure Laminate Core Sidewalls

2. CORE CELL AREA: ≤3 cm<sup>2</sup>

3. WEIGHT: 56.00 lbs

4. ROHS: COMPLIANT



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY



FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING

THIRD ANGLE _ PROJECTION	<b>\$</b>	TITLE	23" x 35", English Lab Table	
ALL DIMS IN	in	DWG NO	54255	SHEET 1 OF 1