
NOTES:

- MATERIAL: STAINLESS STEEL AND
 BLACK ANODIZED ALUMINUM
- 2. MAX TRAVEL DISTANCE: 13mm
- 3. RoHS: EXEMPT



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING

| 81 | | Edmund Optics® |
|----|-----|---|
| | 1.0 | TD 4. / (F) F) F 0.00 F) TUDU 1.01 F 0.41 |

| THIRD ANGLE PROJECTION | | TITLE | 13mm TRAVEL, FINE SCREW, THRU-HOLE BALL BEARING STAGE | |
|------------------------|----|--------|--|-----------------|
| ALL DIMS IN | in | DWG NO | 56356 | SHEET 1 OF 1 |