


NOTES:

- MATERIAL: STAINLESS STEEL AND
 BLACK ANODIZED ALUMINUM
- 2. MAX TRAVEL DISTANCE: .5in
- 3. RoHS: COMPLIANT







SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

FOR INFORMATION ONLY: DO NOT MANUFACTURE PARTS TO THIS DRAWING

	Edmund	Optics®
--	--------	---------

THIRD ANGLE PROJECTION		TITLE	0.5in TRAVEL, ENG MICROMETER, THRU- BALL BEARING STAGE	HOLE
ALL DIMS IN	in	DWG NO	37957	SHEET 1 OF 1