
## NOTES:

- MATERIAL: STAINLESS STEEL AND
   BLACK ANODIZED ALUMINUM
- 2. MAX TRAVEL DISTANCE: 19mm
- 3. RoHS: EXEMPT



Edmund Optics®

∘∥⊓

**-** 22.23

R6.75

DO NOT MANUFACTURE
PARTS TO THIS DRAWING

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE DIMENSIONS ARE FOR REFERENCE ONLY

THIRD ANGLE PROJECTION		TITLE	19mm TRAVEL, FINE SCREW, THRU-HOLE BALL BEARING STAGE	
ALL DIMS IN	in	DWG NO	55285	SHEET 1 OF 1