## NOTES:

- 1. REFLECTIVE MATERIAL: Spectralon®
- 2. NOMINAL REFLECTANCE: N/A
- 3. COLOR: Red, Green, Blue, Yellow
- 4. HOUSING MATERIAL: Delrin
- 5. TEMPERATURE: -80°C 350°C
- 6. RELATIVE HUMIDITY: 5% 95%
- 7. REFLECTANCE STANDARD SHOWN BELOW AS REFERENCE FOR REFLECTANCE STANDARD DIMENSIONS. PART 56079 CONSISTS OF A PREMIUM WOODEN CASE WITH 4 STANDARDS INCLUDED WITHIN.

