



NOTES:

1. REFLECTIVE MATERIAL: Spectralon®
2. NOMINAL REFLECTANCE: 99%, 75%, 50%, 2%
3. COLOR: White and Grayscale
4. HOUSING MATERIAL: Delrin
5. TEMPERATURE: -80 to 350
6. RELATIVE HUMIDITY: 5% - 95%
7. REFLECTANCE STANDARD SHOWN BELOW AS REFERENCE FOR REFLECTANCE STANDARD DIMENSIONS. PART 13513 CONSISTS OF A PREMIUM WOODEN CASE WITH 4 STANDARDS INCLUDED WITHIN.



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE  
DIMENSIONS ARE FOR REFERENCE ONLY

**FOR INFORMATION ONLY:  
DO NOT MANUFACTURE  
PARTS TO THIS DRAWING**

**EO**® **Edmund Optics**®

THIRD ANGLE PROJECTION 	TITLE	Certified Set of 4, Spectralon® White and Grayscale Diffuse Reflectance Standards (99%, 75%, 50%, and 2%)	
	ALL DIMS IN	mm	DWG NO 13513 SHEET 1 OF 1