





NOTES:

1. MATERIAL: BLACK ANODIZED ALUMINUM

2. ROHS COMPLIANT



**FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING**

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE
DIMENSIONS ARE FOR REFERENCE ONLY

EO® Edmund Optics®

THIRD ANGLE PROJECTION		TITLE	50mm, Multiple Filter Holder	
ALL DIMS IN	mm	DWG NO	55007	SHEET 1 OF 1

SIZE OF COMPATIBLE OPTICS (mm)	$\varnothing 50.00$
MAX THICKNESS OF COMPATIBLE OPTICS (mm)	5.0